

Chlorophyll Determination:

Time required: 2day

Procedure

1. Weight 100 mg of leaf tissue in fractions into vial containing 7ml Dimethyl sulphoxide
2. Chlorophyll will extract into the fluid without grinding at 65 C by incubating for various times, depending on the degree of cutinization and thickness of the leaf
3. Transfer the extract liquid to a graduated tube and made up to a total volume of 10 ml with Dimethyl sulphoxide, assay immediately or transferred to vials and stored between 0-4 C until required for analysis
4. Take 3 ml of chlorophyll extract and transfer to cuvette,
5. Measure the optical density (OD) of the extract at the following wavelengths 645 and 663 nm using Dimethyl sulphoxide as a blank after 30 min, 1 hr, incubation.

Calculate total chlorophyll as mg/g of tissue, using the following equations:

$$\text{Chlorophyll A (mg/g)} = 12.7 (\text{OD}_{663}) - 2.69(\text{OD}_{645}) \times (V/(1000 \times \text{wt}))$$

$$\text{Chlorophyll B (mg/g)} = 22.9 (\text{OD}_{645}) - 4.68(\text{OD}_{663}) \times (V/(1000 \times \text{wt}))$$

$$\text{Total Chlorophyll (mg/g)} = 20.2 (\text{OD}_{645}) + 8.02(\text{OD}_{663}) \times (V/(1000 \times \text{wt}))$$

Where;

OD: optical density at certain wave length (645 or 663 nm)

v: final volume (10 ml)

wt: weight of sample (100 mg)