

# Measurement of pH and Titratable Acidity

**Time required 1day**

## **Procedure**

- A. Obtain at least 50 mls of clear juice by one of the following methods:
1. Cut fruit, press with a hand press, and filter through cheesecloth, or
  2. Cut fruit into a blender, homogenize, centrifuge slurry, and pour off clear liquid for analysis.
- \*\* Sugar levels often vary within the fruit, being higher at the stem-end and lower at the calyx-end. For this reason, it is important to use longitudinal slices of fruit (from end to end) when sampling.
- B. Make sure samples are at room temperature before taking measurements.
- C. Measure the pH of the samples with a pH meter and record the value.
- D. For each sample, weigh out 6 grams of juice into a 100 ml beaker.
- E. To each sample, add 50 mls of water.
- F. Titrate each sample with 0.1 N NaOH to an end point of 8.2 (measured with the pH meter or phenolphthalein indicator) and record the milliliters (mls) of NaOH used.

**G. Calculate the titratable acidity using the following formula:**

$$\% \text{ acid} = \frac{[\text{mls NaOH used}] \times [0.1 \text{ N NaOH}] \times [\text{milliequivalent factor}] \times [100]}{\text{grams of sample}}$$